

DKT. 13131

17608

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

Applicant: Susan L. Weston et al.

Art Unit:

1655

Serial No.:

09/228,639

Examiner:

Enewold, J.

Date Filed:

12 January 1999

Docket No.:

13131

For:

SEQUENCES

September 10, 2002

Hon. Commissioner of Patents and Trademarks Washington, D.C. 20231

SUPPLEMENTAL DECLARATION

Sir:

- I, Dr. Gary Brown, am skilled in the art of Molecular Biology. I am a scientific
 consultant in biotechnology and am being compensated for making this declaration.
- I make this declaration based upon my training, knowledge, education, and experience as a molecular biologist, my review of the application and the history of the prosecution of this application as reflected in the file maintained by the attorneys prosecuting this application, and my review of prior art cited by the Examiner during prosecution of this application, and my review of inventor laboratory monthly project reports (laboratory reports) and the CF 12 Test Prototype Record Sheets (record sheets) attached as Exhibit D.
- 3. In the March 12, 2002 Office Action, the Examiner states on page 3 that "[T]he exhibits of the laboratory reports however do not appear to identify primers by either

Susan L. Weston et al.

Serial No.: Filed: page 2 09/228,639 12 January 1999

sequence or by SEQ ID NO: as relied upon in the instant specification. The exhibits are general to primers designated by mutational name and do not provide any information as to the structure of the primers which fail. Therefore, it is unclear whether these primers which failed are the primers claimed or whether they are within the scope of the claims or whether they constitute unexpected results. Further, the declaration asserts that primer length affected the results of the assay. As explained above, while the declaration asserts primers of different lengths were used, it is not clear from the manual or from the declaration which of these lengths was ideal and whether this result is within the scope of the claims."

- 4. In this declaration, I provide the record sheets that I reviewed, which identifies, among other things, some of the primer sequences used by name, changes made to primer sequences, nucleotide lengths and concentrations.
- Based upon my review of, among other things, the laboratory reports (previously submitted as Exhibits B and C) of Susan Weston (Susie Weston) and the record sheets, as one skilled in the art, I conclude that the laboratory reports and record sheets detail the unexpected results of the invention. The reports and record sheets record difficulties, the multiple failures, and the surprising results of the multiple experiments performed by the inventor in the present application. These experiments show unexpected results after modifications of numerous parameters of the experimental system.
- 6. Inventor Susie Weston, in her laboratory report from Lab books NBY 5882 and NBY.

Susan L. Weston et al.

Serial No.: Filed:

09/228,639 12 January 1999

page 3

5922, and record sheets shows that certain primers were prone to non-specific priming which makes such primers unacceptable for use in this assay. A number of unacceptable primer dimer results were observed, and a reversed direction primer was selected and introduced as a standard, to avoid primer dimers (Exhibit B, page 1). In addition, zinc concentrations and temperature were varied. Some of the unacceptable results the inventor documented were weak signal, very weak signal, primer dimer, and nonspecific bands (Exhibit B, page 2).

In further experiments, inventor Susie Weston designed primers (see for example, 7. record sheets documenting changes to prototype 2 and prototype 3, involving primers: KRT 16, R334WM, R1162XM and R1162XC), varied concentrations of primers, tried zinc-free preparations, and varied the annealing times. Results of these modifications that the inventor recorded included unacceptably faint signals and smeared products, absent product, and nonspecific products (Exhibit B, page 3). For example, the inventor found that primer R553XM (see prototype 1 of the record sheets) did not perform as expected and changed the 3' sequence of the primer from CCA to TGA. The changed sequence was found to work well in conjunction with the other primers and was introduced as a standard in the set and is now designated as SEQID No. 24. Further the inventor found the primers KRT 16 reverse, R334WC forward and R1162XM reverse (see prototype 2 of the record sheets) did not function and were changed dramatically. For example, in Exhibit B the inventor states "[t]aking into account that the changes in mixes 2 and 3 had no adverse effect and were potentially beneficial the new 621+1/R117HC, 0.5μ M R334WC, R1162XM and R1162WC and 2μ M621+1M and R553XM were introduced as standard." These primers introduced

Susan L. Weston et al.

Serial No.: Filed: 09/228,639 12 January 1999

page 4

as standard are named: R334WC (SEQ ID No:22), R1162XM (SEQ ID NO: 26) and R553XM (SEO ID No: 24). Primer 621+1/R117HC was changed when it failed to work as expected in the final stages of prototype testing (see last prototype CF(12)m Test: Prototype Record Sheet B Tube [P8] of record sheets). This new 621+1/R117HC primer is now (SEQ IN No: 20). Interestingly, the inventor later found the new R1162WC surprisingly did not function well (see Exhibit C, first page, last paragraph) and reinstated primer R1162XC as standard now SEQ ID NO: 27. Accordingly, one of ordinary skill in the art can easily see from the sheer number of modifications to the assays (primer sequence and chemical modifications) that the inventor did not simply optimize a functioning assay. Rather, the inventor had to discover an assay and all of its components by changing primer sequences, primer concentrations and chemical components of the assay in a way that was not routine to achieve the claimed invention. Thus, as can be seen in the specification, notebook pages and record sheets, I conclude that this inventor designed multiple variations of the experiment expecting to overcome problems of nonspecificity, weak products, and primer dimer contaminants, she observed unexpected results and had to discover and assay that worked by numerous trials, drawing on further variations and further experiments that were not routine to achieve the claimed invention.

8. Inventor Susie Weston, in her laboratory report from Lab books NBY 5922 and NBY 5935 (Exhibit C) concludes that the results she obtained from further experiments were unexpected and were not routine. Following experiments that attempted to eliminate nonspecific products and increase band intensities—problems from the previous months' results—some resulting signals were again too weak, and so she

Susan L. Weston et al.

Serial No.: Filed: page 5 09/228,639 12 January 1999

varied the concentrations of primers. She also varied primer length, using 29mers, 26mers, 24mers, and also temperature changes in the protocol as well (Exhibit C, page 1). I particularly note the following observation of the R1162XC primer (see primer sequence on prototype 3 record sheet) by her: "This surprisingly had the reverse effect -all ...DNAs were specific with the original primer and non-specific with the new primer" (Exhibit C, page 1, emphasis added). She went on to vary concentration of primers, lengths of primer, dropped certain primers, varied annealing temperatures, varied Taq enzyme units, and varied DNA concentrations (Exhibit C, page 2). Results she described were weak, very weak, nonspecific—i.e. unacceptable.

- 9. Inventor Susie Weston, in her laboratory report from Lab books NBY 5935 and NBY 6107, shows her experimental variation of dNTP reagent concentration, and design of new primers to attempt to overcome mispriming difficulties. Primer length and identity of primers were varied. Her results included nonspecific bands, weak 'diagnostics, and extremely weak non-specific bands (Exhibit C, page 3). A result noted by the inventor states "[i]t was thought it might come up 621+1 as only the 3' base would be mismatched but much to our relief it didn't." (Exhibit C, page 3, paragraph 3). Nonspecific signals and weak signals continued to be problematic, which shows that even more experimentation was needed. Accordingly, I conclude that the inventors modifications sometimes resulted in improved assays other times resulted in new or continued problems, and that these results to obtained the presently disclosed primers were not routine.
- 10. At page 3 of the March 12, 2002 Office Action, the Examiner states: "[als noted in In

Susan L. Weston et al.

Serial No.: Filed: page 6 09/228,639 12 January 1999

re Aller, 105 USPQ 233 at 235, more particularly, where the general conditions of a claim are disclosed in the prior art, it is not inventive to discover the optimum or workable ranges by routine experimentation." The Examiner further states that routine optimization is not considered inventive and no evidence has been presented that primer selection performed was other than routine, that the products resulted from optimization have any unexpected properties, or the results should be considered unexpected in any way as compared to the closest prior art.

- 11. I respectfully disagree with the Examiner. None of the cited references relied upon by the Examiner disclose the primer sets of the claimed invention that allow optimization and production of meaningful results in the ARMS assay. In my opinion, the presently claimed primer sets were not known prior to the instant invention, and it is not routine optimization to select these primer sets out of the numerous primers that did not work and get them to function in the ARMS assay. As shown by the specification, the exhibits submitted and the declarations, Applicants did not simply optimize a known set of primers or ARMS parameters. Applicants discovered how to design primers of the instant invention and establish conditions to allow those primers to function together simultaneously in primer sets to yield reproducible and valid results. Accordingly, this was not routine experimentation.
- 12. I hereby declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true, and further that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001

Susan L. Weston et al.

Inventor: Serial No.: Filed:

09/228,639

page 7

12 January 1999

of Title 18 of the United States Code and that such willful false statements may jeopardize the validity of the application or any patent issued thereon.

Dr. Gary L. Brown

)ate

ZENECA Diagnostics

F12 Test: Prototype Record Sheet

Tube Prototype	Operator SUSIE	Date	Standard ARMS/dNTP?	ues.
			•	

	•
Modifications since last prototype	

Primer	F/R	Oligo ID	Mismatch	Length	Conen (µM)	Sequence	
APOB	F	29.500	C-A-2	29	0.25	CTCTGGGAGCACAGTACGAAAACCAC CT	
Apo B	R	004895	C-A-2	30	0.25	CAGCCAAAACTTTTACAGGGATGGAGAA	CG
KRT16	F	144990	C-T-2	28	0.1	CTGAGCACCTTTCCTTC.TTTCAGCCACC	
KRT16	R	M4991	A-C-2	28	0.1	CATGGAGCTGGAGGAGGTGAACTGGCAG	<u>†</u>
AF N	R	0017-94-	C-T-2	30	(GTATCTATATTCATCATAGGAAACACC AC	A
△FC-	F	001894	1	30		GACTICACTICTAATGATGATTATGGGAG	rΑ
62.HIM	R	2671	G-T-2	30	1	TGCCATGGGGCCTGTGCAAGGAAGTATTTG.	A.
R117HM	R	245411	C-C-2	30	Ľ	AGCCTATGCCTAGATAAATCGCGATAGAC	T
621/RIJAH	F	2678P	_	30	l	TCACATATGGTATGACCCTCTATATAA AC	T
R334WM	R	245402	A-G-2	30	1	GCAGAATGAGATGGTGGTGAATATTTTCA	A
R334WC	F	245404		30		GAAGGCAGCCTATGTGAGATACTTCAA TAG	G
G5510M	R	2020.94	T-T-2	ఎక	1	GCTAAAGAAATTCTTGCTCGTTGTT	
2553XM	R	001694	C-A-3	60	i	GACTGACTGACTGACTCTGACTGACT	7H
				nhtail		TCACCTTGCTMAGAHATTCTTGCCCA	4
CYSKINDERS SK	E	0043.95	_	28		TAAAATTTCAGCAATGTTGTTTTT GA	+ CC
R1162XM	R	245409	A-G-2	30		TGTTGGCATGTCAATGAACTTAAAGAC	TA.
RIIBZXC	۳	245405		3O	(AATTATTTCCTGTTAGTTCATTGAAAAGC	C

Process Number	Quantity	Destination	Comments
4145-06	60	Mike W.	sent
			94°C-1,58°C-1,72°C-1, x35, 2UTag/test.
4145.06	60	RMCH	sent
			94°C-1'58°C-1'72°C-1'x35 QUTag/test.
·			

Other Information:

cf12prot.sam



Tube B Prototype 2 Operator SUSIC Date Standard ARMS/dNTP? 465 Modifications since last prototype [APOB], [KRT], R553XM mismatch

Primer	F/R	Oligo ID	Mismatch	Length	Сопсп (µМ)	Sequence	
APO B	F	0047-95	C-A-2	29	0.1	CTCTGGGAGCACAGTACCTATTA AACCAC	C.T
Apo B	R	0048.95	C-A-2	30	0.1	CAGCCAAAACTTTTACAGGGTATGGAGA	ACG
KRTIL	F	M4-990	C-T-2	28	0.25	CTGAGCACCTTICCTTCTTCAGCCACC	
KRTIB	R	M4991	A-C-2	28	0.25	CATGGAGCTGGAGGAGGTGAACTGGCA	Gr
∆FN	R.	0017-94	C-T-2	<i>30</i>	1	GTATCTATATTCATCATAGG A MACACCAC	ţ A
AFC	F	0018.94		30		GACTICACTICTAATGATGATTATGGGAC	711
621+1M	R	2671	G-T-2	30	1	TGCCATGGGGCCTGTGCAAGGAAGTAT	TGA
RITHIM	R	245411	C-C-2	30		AGCCTATGCCTAGATAAATCGCGATAG	ACT
GZI/RIAH	F	2678P	٠	30	1	TCACATATGGTATGACCCTCTATATAAA	CT
2334WM	R	24-5402	A-G-2	30)	GCAGAATGAGATGGTGGTGAATATTT	TCAA
R334WC.	F	24-5404		30	1	GAAGGCAGCCTA TGTGAGATACTICA	ATAG
GSSIDM	R	0020-94	ナナ2	25	1	GCTAAA GAAATTCTTGCTCGTTGTT	
RSSEXM	R	ಎಚಿತ	G-G-2	60	1	GACTGACTGACTGACTCTGACTG	ACT
				nhicui	· .	TATTICACCTTGCTAAAGAAATTCTTGC	TGA
GSSI DIESSA C	F	0043-95		28	1	TAAAATTICAGCAATGTTGTTTTGAC	C
RIIBZXM	R	245409	A-G-2	30)	TGTTGGCATGTCAATGAACTTAAAGAC	TAA
RII62XC		245405		30	1	AATTATTTCCTGTTAGTTCATTGAAAA	GCC.

Mix	es	sent	out:

Process Number	Quantity	Destination	Comments
2175-03	?	MIKE W.	sent
			94°C-11, 59°C-1, 72°C-3050C x35 2UTaq / tsl
4185-03	144	MIKEW.	sent
			94°C-1, 58°C-1, 72°C-1, 35 2UTag/rest

Other Information:

4185-03 contained a repeat primer syntheses: 0017-94 replaced by 0058-95

cois-94 replaced by ocs9-95. cf12prot.sam



Tube B Prototype 3 Operator SUSIC Date Standard ARMS/dNTP? UES

Modifications since last prototype 10+5 & CP2 Was April)

Primer	F/R	Oligo ID	Mismatch	Length	Concn (µM)	Sequence
APOB	F	0159-95	C-A-2	29	0.1	CTC, TGG, GAG, CAC, AGT, ACG, AAA, AAC, CAC CT
APOB	R	0160-95	C-A-2	<i>3</i> Q	0.1	CAG CCA AAA CTT TTA CAG GGA TGG AGA ACG
KRT16	F	0161-95	C-T-2	28	0.25	CTG AGC ACCTTTCCTTCTTTC AGC CAC C
KRT16	R	0162.95	A-C-2	27	0.25	TTC ATG GAG CTG GAG GAG GTG AAC TAG
DF N	R	0176-45	C-T-2	30)	1	GTATCTATATTCATCATAGGAAACACC ACA
DF C	F	0177-95		30		GACTICACTICTAATGATGATTATGGG AGA
621+1M	R	0178-95	G-T-2	30	2	TGCCAT GGG GCCTGTGCA AGG AAGTAT TGA
R117HM	R	0179.95	C-C-2	30	1	AGCCTATGCCTAGATHAATCGCGATAG ACT
621/R117H	F	0196-95	C-7-6 T-C-26 T-3-27	34	1	TCACATAGTGTATGACCCTCTATATACACTCATT
R334WM	F	P1337	G-G-2	30	0.5	CLTATGCACTAATCAAAGGAATCATCCTGT
R334WC	R	0192-95		30	0.5	TTTGTTTATTGCTCCAAGAGAGTCATA CCA
GSSIDM	R	0190-95	T-T-2	25		GCT AAAGAA ATTCTTGCTCGTTGTT
RSS3XM	R	0191-95	G-G-2	60	ュ	GACTGACTGACTGACTGACTGACTG
				nhtail		TATTCACCTTGCTAAAGAAAT TCTTGCTG
45210/1880	F	0193.95	G-G-20 G-A-21	28	1	TAAAATTGGAGCAATGTTGTTTT GACC
R1162XM		0157-95		30	0.5	TATTITITITICAGATGCGATCTGTGA GTT
R1162XC	R	P5732		29	0.5	TTTTGCTGTGAGATCTTTGACAGTCATTT

Mixes sent out:

Process Number	Quantity	Destination	Comments
4345-03	50	RMCHI.	sent Noresults after 6
			weeks so sent next prototype
			, , , , , , , , , , , , , , , , , , , ,
	-		
	ł	ł	·

Other Information:

cycling: 94°C-1,58°C-2,42°C-1' 24 Tag | reaction. Recommended 50ng DNA/reaction.

cf12prot.sam



Tube Prototype 4 Operator SUSIC Date Standard ARMS/dNTP? Mes Shorter Apo B primers. dec. EAFJ. dec. ERIIBZX]

Modifications since last prototype olec. [R334W] inc. [GSSID]RSSBXC] inc [621+1/R117+11C]

inc. EKRTJ

Primer	F/R	Oligo ID	Mismatch	Length	Concn (µM)	Sequence :	
Apo B	F	0197-95	C-A-2	24	0.1	GGAGCACAGTACGAA AAACCACCT	
APCB	R	0198.95	C-H-2	26	0-1	CAA AACTITTACAGGGATGGAGAA C	G
KRT16	F	0161-95	C-T-2	28	0.3	CTGAGCACCTTTCCTTCTTTCAGCCAC	\subset
KRTIB	R	0162-95	A-C-2	27	0.3	TTCATGGAGCTGGAGGAGGTGAACTA	1G-
AFN	R	0176-95	C-7-2	30	0.5	GTATCTATATTCATCATAGGAAACACCA	1CA
AFC	F	0177-95	-	30	0.5	GACTICACTICTAATGATGATTATGGGA	4 <i>G</i> -A
621+1M	R	0178-95	G-T-2	30	a a	TGC CATGGGGCCTGTGCAAGGAAGT	ATTG
RIITHM	R	0179.95	C-C-2	30	1	AGCCTATGCCTAGATAAATCGCCATAG	ACT
621/21171-1 C	F	0196-95	C-T-8 J-5-24	34	1.5	TCACATAGTGTA TGACCC.TCTATATACACTC	ATT
R334 WM	F	P1337	G-G-2	30	0.25	CCTATGCACTAATCAAAGGAATCATCC	167
R334WC	R	092-95		30	0.25	TITGITTATTGCTCCAAGAGAGAGTCATAC	ĽΑ
G551D14	R	0190.95	T-T-2	25	1	GCTAAAGAAATTCTTGCTCGTTGTT	1
RSS3XM	R	0191-95	G-G-2	60	<u>ي</u>	GACTGACTGACTGACTGACTGACTGACT	77A7
				nhtail		TCACCTTGCTAAAGAAATTCTTGCTGA	
75510/12553X C	F	0193.95	G-G-20 G-A-21	૨૬	1.5	TAMATTGGAGCAATGTTGTTTTTGACC	-
P116ZXM	1	0157-95		<i>3</i> 0	0.25	TATTTTTATTTC AGATGCGATCTGTGAG	77
R1162XC.	R	P5732		29	0.25	TTTTGCTGTGAGATCTTTGACAGTCAT	TT

Mixes sent out:

Process Number	Quantity	Destination	Comments
1395-01	2004	MIKKI. D.	
1395-01	22	Nikki D.	sent to evaluate alkali
			+ phenol extraction cfcf blocals
139501	50	RTCH	sent
1395-01	20	MIKKI D.	sent .

Other Information:

Cycling 94°C-1; 58°C-2! 72°C-1' 247aq/reaction Recommended 50ng DNA/reaction.

cf12protsam



Tube B Prototype 5 Operator SUSIE Date Standard ARMS/dNTP? 4eS

Modifications since last prototype ODC CONTYCLS, CleC. EthpoBJ.

Primer	F/R	Oligo ID	Mismatch	Length	Concn (µM)	Sequence
APOB	F	0197.95	C-A-2	24	0.075	GGAGCACAGTAC GAAAAACCACCT
APCE		0198-95		26	0.075	CAAAACTTTTACAGGGATGGAGAACC
ODC		0235-95		30	0.3	AGAGGATTATCTATGCAAATCCTTGTA
ODC	R	0236-95	A-C-2	30	0.3	TCAACTTCACTATCAAAAGTCATCATC
AFN	R	0176.95	C-T-2	30	0.5	GTATCTATATTCATCATAGGAAACACC
ΔFC	F	0177-45		30	0.5	GACTICACTICTAATGATGATTATGGG
621+1M	R	0178-95	G-T-2	30	2	TGCCATGGGGCCTGTGCAAGGAAGT
RIFHM		0179-95		30	}	AGCCTATGCCTAGATAAATCGCGA"
621 BITH	F	096-95	C-1-36	34	1.5	TCACATAGTGTATGACCCTCTATATACAC
R334WM	F	P1337	G-G-2	<i>3</i> 0	0.25	CCTATGCACTAATCAAAGGAATCATCCTC
2334WC	R	0192-95	_	30	0.25	TITGTTTATTGCTCCAAGAGAGTCATA
GSSIDM	R	0190-95	T-T-2	25		GCTAAAGAAATTCTTGCTCGTTGTT
PS53XM	R	019195	G-G-2	60	2	GACTGACTGACTGACTGACTGACTGA
				nHail		TATTCACCTTGCTAAAGAAATTCTTGC
CLERIDISIZIX	F	019345	G-G-20 G-A-21	28	1-5	THAHATTGGAGCAATGTTGTTTTG,
E1162XM	F	0157.95	G-T-2	30	0.25	TATTTTATTTCAGATGCGATCTGTG
P1162XC	R	P5732		29	0.25	TTTTGCTGTGAGATCTTTGACAGTCAT

Process Number	Quantity	Destination	Comments
4465-02	53	Nikki D.	Box 4 for 50 mouthwashes, sent
3435-04-	120	KGcentre	tock
3495-04	55	MIKKID.	for 50 mouthwashes sent
		71	using fresh control primer syntheses.
			·

Other Information:

Cycling: 94°C-1, 58°C-2, 72°C-1 ×35

24 Tag/reaction

of 2010 sam 50ng DNA/reachon or alk prop



Tube B Prototype O Operator SUSIC Date Standard ARMS/dNTP? UCS

Reven Shorker Apoils (at the conch!)

Modifications since last prototype dec. CRITHING CRITEXING CRITEXING CRITEXING CRITEXING

CRITERIAN

CR

	Primer	F/R	Oligo ID	Mismatch	Length	Concn (µM)	Sequence
45	Apo B	F	0304-45	C-A-2	ನ್ನು	0.1	AGC ACA GTA CGA AAA ACCACCT
ĺ	Apos	R	0305-95	C-A-2	ゔゔ	0.1	ACT TTTACA GGG ATG GAG AACG
	Ω Ω	Ш	0302-95	C-T-2	30	0.3	AGAGIGIATTATC TATGCAAATCCTTGTA ACC
	X	R	0303-95	A-C-2	30	0.3	TCAACTTCACTATCAAAAGTCATCATCT AG
	△FN	Ù	0176-95	C-T-2	30	0.5	GTATCTATATTCATCATAGGAAACACC ACA
	∆FC.	Į.	0177-95		SC.	0.5	GACTTC ACTTCTAATGATGATTATGGG AGA
	621M	R	0175.95	G-T-2	30	.⊋	TGCCATGGGGCCTGTGCAAGGAAGTATTGF
	21711	R	0+79-95	<u>C-C-2</u>	30	Ċ Ċ	AGCCTATGCCTAGATAAATCGCGATAG ACT
	621/P1741	F	0196-95	ر در در در این در در این	34	<u>س</u> ن	TCACATAGTGTATGACCCTCTATATACACTCAT T
	R334WM	Ţ	171337	G-G-2.	30	<u>ر</u> ن	CCTATGCACTAATCAAAGGAATCATCCTGT
	P334WC	R	0192-95		3O	0.2	TITGTTTATTGCTCCAAGAGAGTCATACCA
	GSSIDM	R	0190.95	下T-2	んい		GCTAAAGAAATTCTTGCTCGTTGTT
	RSSISKM	R	0301-45	G-G-2	60	3	GACTGACTGACTGACTGACTGACTGACTTAT
			·	·——	nhtail		TCACCITGCTAAAGAAATTCTTGCTGA
	USSIDJESSKY C	F	0193-95	G-G-20 G-A-21	28	1.5	TAAAATTGGAGC AATGTTGTTTTTGACC
	RIGZKM		0157.95		30	0.2	TATTTTTTTCAGATGCGATCTGTGAGTT
	RII62XC	R	P5732		ত্ৰপ	0.2.	TTTTGCTGTGAGATCTTTGACAGTCATTT

Mixes sent out:

Process Number	Quantity	Destination	Comments
	-		•

Other Information:

now 2-5 litting/reaction.

cf12prot.sam

B Tube [P7]

1x ARMS, 100µM dNTPs

Modifications since last prototype:

new ApoB primers, each at 0.05 MM

Primer	F/R	Oligo ID	Mismatch	Length	Concn (µM)	Sequence
Аро В	F	S5890	CA -2	23	0.05	GAG CAC AGT ACG AAA AAC CAC CT
Аро В	R	S5891	CA -2	24	0.05	AAA CTT TTA CAG GGA TGG AGA ACG
ODC	F	0302-95	CT-2	30	0.30	AGA GGA TTA TCT ATG CAA ATC CTT GTA ACC
ODC	R	0303-95	AC-2	30	0.30	TCA ACT TCA CTA TCA AAA GTC ATC ATC TAG
DF N	R	0176-95	CT-2	30	0.50	GTA TCT ATA TTC ATC ATA GGA AAC ACC ACA
DF C	F	0177-95	-	-30	0.50	GAC TTC ACT TCT AAT GAT GAT TAT GGG AGA
621 M	R	0178-95	GT-2	30	2.00	TGC CAT GGG GCC TGT GCA AGG AAG TAT TGA
R117H M	R	0179-95	CC -2	30	0.50	AGC CTA TGC CTA GAT AAA TCG CGA TAG ACT
621/R117H C	F	0196-95	CT -8 TC -26 GA -27	34	1.50	TCA CAT AGT GTA TGA CCC TCT ATA TAC ACT CAT T
R334W M	F	P1337	GG -2	30	0.20	CCT ATG CAC TAA TCA AAG GAA TCA TCC TGT
R334W C	R	0192-95	-	30	0.20	TTT GTT TAT TGC TCC AAG AGA GTC ATA CCA
G551D M	R	0190-95	TT-2	25	1.00	GCT AAA GAA ATT CTT GCT CGT TGT T
R553X M	R	0301-95	GG -2	60	2.00	GAC TGA CTG ACT GAC TGA CTC TGA CTG ACT TAT TCA CCT TGC TAA AGA AAT TCT TGC TGA
G551D/R553X C	F	0193-95	GG -20 GA -21	28	1.50	TAA AAT TGG AGC AAT GTT GTT TTT GAC C
R1162X M	F	0300-95	GT -2	30	0.20	TAT TIT TAT TIC AGA TGC GAT CIG TGA GTT
R1162X C	R	P5732	-	29	0.20	TTT TGC TGT GAG ATC TTT GAC AGT CAT TT

Operator susie weston Date

ZENECA Diagnostics

B Tube [P8]

1x ARMS, 100µM dNTPs

Modifications since last prototype:

621+1/R117HC 3 bases larger at 5'end

Primer	F/R	- 1.5-1.5	Mismatch	Length	Concr (µM)	Sequence
Аро В	F	S5890	CA -2	23	0.05	GAG CAC AGT ACG AAA AAC CAC CT
Аро В	R	S5891	CA -2	24	0.05	AAA CTT TTA CAG GGA TGG AGA ACG
ODC	F	0302-95	CT-2	30	0.30	AGA GGA TTA TCT ATG CAA ATC CTT GTA ACC
ODC	R	0303-95	AC -2	30	0.30	TCA ACT TCA CTA TCA AAA GTC ATC ATC TAG
DF N	R	0176-95	CT-2	30	0.50	GTA TOT ATA TTO ATO ATA GGA AAC ACC ACA
DF C	F	0177-95	-	30	0.50	GAC TTC ACT TCT AAT GAT GAT TAT GGG AGA
621 M	R	R029-96	GT-2	30	2.00	TGC CAT GGG GCC TGT GCA AGG AAG TAT TGA
R117H M	R	0179-95	CC -2	30	0.50	AGC CTA TGC CTA GAT AAA TCG CGA TAG ACT
621/R117H C	F	S7587	CT -8 TC -26 GA -27	37	1.50	GTT TCA CAT AGT GTA TGA CCC TCT ATA TAC ACT CAT T
R334W M	F	P1337	GG-2	30	0.20	CCT ATG CAC TAA TCA AAG GAA TCA TCC TGT
R334W C	R	0192-95	-	30	0.20	TIT GTT TAT TGC TCC AAG AGA GTC ATA CCA
G551D M	R	0190-95	Π-2	25	1.00	GCT AAA GAA ATT CTT GCT CGT TGT T
R553X M	R	0301-95	GG -2	60	2.00	GAC TGA CTG ACT GAC TGA CTC TGA CTG ACT TAT TCA CCT TGC TAA AGA AAT TCT TGC TGA
G551D/R553X C	F	0193-95	GG -20 GA -21	28	1.50	TAA AAT TGG AGC AAT GTT GTT TTT GAC C
R1162X M	F	0300-95	GT-2	30	0.20	TAT TIT TAT TIC AGA TGC GAT CTG TGA GTT
R1162X C	R	P5732	-	29	0.20	TTT TGC TGT GAG ATC TTT GAC AGT CAT TT

Operator 848ie Westou

Date

Note: Oswel synthesis S7587 gave primer-dimer. It was replaced by Nikki's synthesis R072-96 which was fine.

ZENECA Diagnostics